


# MR 3000 Manual Rotary Microtome




The sample thickness can be set easily by using the precision setting knob




Precise orientation mechanism makes it easy to adjust the specimen to the right position

Coarse feed wheel

Handwheel can be locked in any position



Lateral blade holder adjustment. The entire width of blades can be used without releasing the blade



Universal fast release cassette clamp, suitable for all type of embedding cassettes or rings



Large sections waste tray

Want More Information ? [www.histo-online.com](http://www.histo-online.com)



## Technical Specifications

- Section thickness range: 0 ~ 60µm
- Setting value: 0 ~ 20µm increment 1µm  
20 ~ 60µm increment 5µm
- Minimum Setting value: 1µm
- Precision Error: ±5%
- Maximum Specimen Size: 60x50mm
- Trimming Setting Mode: 15µm or 35µm
- Specimen Retraction Range: 90µm
- Specimen retraction range of 90µm can be shut off
- Object Feed: 30mm
- Vertical Stroke: 70mm
- Specimen Orientation System: 8 X/Y
  
- Dimension: 570mmx500mmx300mm
- Net Weight: 33kg




Distributed by:



Via Brembo, 27 - Milano Italy Tel: +39 0255230061 Fax: +39 0255213764 Email: [histoline@histoline.com](mailto:histoline@histoline.com) URL: [www.histoline.com](http://www.histoline.com)

NB: Specifications are subject to change without any prior notice



# MR 3000

Manual Rotary Microtome  
with accurate trimming function  
Safe, Comfortable, Reliable

