


ATS 200




Automatic Slide Stainer

for Multiple Histological & Cytological Routine Stains



ATS 200

ATS 200 is specially designed for simultaneous process of different staining protocols, allows continuous unload & reload without interrupting staining cycle. All functions, menu-driven programming, incubation times and sequence of use of reagent stations are freely programmable and controlled by the touch screen panel. Ergonomic design, integrated oven for slides drying, digital electronic temperature control, movement parts are in bearing linear guide for high precision on axis.



Technical Specifications

Number of Stations:	26
Number of reagents Stations:	18
Number of wash Stations:	3
Number of waiting Stations:	3
Number of programs:	5
Loading capacity per rack:	30 slides, 200 staining per hour
Capacity of single tray:	450 ml
Incubation time setting:	From 1 min. to 99 h 59 min.
Temperature range:	From ambient to 75°C increments 1°C
UPS:	Optional
Supply Voltage:	220 V / 50 Hz 110 V / 60 Hz
Power draw:	300 V/A
Dimension:	
Width: 1200 mm	Depth: 500 mm
Height: 600 mm	Weight: 110Kg.